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
The image shows a 'New Document Wizard' dialog box with a blue title bar and a close button. The main area is titled 'Define Document' and contains instructions: 'Select the assay, container, and template for the document, and enter the operator name and comments.' Below this are several input fields: 'Assay' (dropdown menu with 'Background' selected), 'Container' (dropdown menu with '96-Well Clear' selected), 'Template' (dropdown menu with 'Blank Document' selected and a 'Browse ...' button next to it), 'Run Mode' (dropdown menu with 'Standard 7500' selected), 'Operator' (text field with 'liangy2' entered), and 'Comments' (text area with 'SDS v1.3.1' entered). At the bottom is a 'Plate Name' field with 'Platel' entered. The bottom of the dialog has four buttons: '< 上一步(B)' (disabled), '下一步(N) >' (active), '完成' (disabled), and '取消' (disabled).

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


The image shows a 'New Document Wizard' dialog box with a blue title bar and a red close button. The main area is titled 'Define Document' and contains instructions: 'Select the assay, container, and template for the document, and enter the operator name and comments.' Below this are several input fields: 'Assay' (Pure Spectra), 'Container' (96-Well Clear), and 'Template' (Blank Document), each with a dropdown arrow. A 'Browse ...' button is next to the Template field. The 'Operator' field contains 'Administrator'. The 'Comments' field contains 'This is an example of a plate document created for a Pure Dye calibration.' and has scroll arrows. The 'Default Plate Name' field contains 'Plate 3'. At the bottom are four buttons: '< Back', 'Next >', 'Finish', and 'Cancel'.

New Document Wizard

Define Document
Select the assay, container, and template for the document, and enter the operator name and comments.

Assay: Pure Spectra

Container: 96-Well Clear

Template: Blank Document

Browse ...

Operator: Administrator

Comments: This is an example of a plate document created for a Pure Dye calibration.

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
The screenshot displays the 7500 Software v2.0.3 interface. On the left sidebar, the 'ROI' status is highlighted in a yellow box, showing 'Status: Current' in green text, 'Last Run: 8/21/09 3:45 PM', and 'Exp. Date: 2/21/10 3:45 PM'. Below this, the 'Background' and 'Optical' sections both show 'Status: Not Calibrated' in red text. The main panel features a 'ROI Calibration' section with instructions and two buttons: 'Start Calibration' and 'Start Manual Calibration'. Below this is a 'ROI Calibration Data' table with five rows of data.

Filter	Status	Last Run	Expiration Date
1	Current	8/21/09 3:49 PM	2/21/10 3:49 PM
2	Current	8/21/09 3:48 PM	2/21/10 3:48 PM
3	Current	8/21/09 3:47 PM	2/21/10 3:47 PM
4	Current	8/21/09 3:46 PM	2/21/10 3:46 PM
5	Current	8/21/09 3:45 PM	2/21/10 3:45 PM

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The screenshot displays the 'Dye Calibration' software window. On the left is a vertical sidebar with three buttons: 'Setup' (highlighted in yellow), 'Run', and 'Analysis'. The main area is titled 'System Dye Calibration Setup' and includes a 'Setup Help' link. It contains the following instructions:

Instructions: Locate the required materials, prepare plate 1, load plate 1 into the instrument, then click "Next."

1. Locate required materials:

- Spectral calibration plates from the spectral calibration kit: plate 1 (FAM™ dye, ROX™ dye, Instrument Normalization) and plate 2 (VIC® dye, JOET™ dye, SYBR® Green I dy...
- Powder-free gloves
- Safety goggles
- Centrifuge with plate adapter

2. Prepare plate 1:

[View Instructions](#)

3. Load plate 1 into the instrument:

- a. Open the instrument drawer.
- b. Place the plate in the sample block with the A1 position at the back-left corner of the sample block.
- c. Close the instrument drawer carefully.

d. Select the checkbox to confirm plate 1 is loaded into the instrument:

☒ Plate 1 (FAM dye, ROX dye, and Instrument Normalization) is loaded into the instrument.

At the bottom right, there are four buttons: 'Cancel', 'Previous', 'Next', and 'Finish'.

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